

## 开发指南

### 1 环境配置

#### 1.1 添加framework

#### 1.2 配置蓝牙权限

### 2 使用

#### 2.1 设备支持的服务UUID

#### 2.2 导入framwrok库

#### 2.3 扫描戒指

#### 2.4 取消扫描设备

#### 2.5 连接戒指

#### 2.6 断开连接

### 3. 戒指支持的指令

#### 3.1 设置手环时间

#### 3.2 读取戒指电量

#### 3.3 绑定戒指亮灯

#### 3.4 设置戒指时间进制/用户个人信息

#### 3.5 获取戒指时间进制/用户个人信息

#### 3.6 获取戒指固件的版本号

#### 3.7 获取当前计步信息

#### 3.8 获取某天总的统计数据

#### 3.9 获取某天的详细运动数据

#### 3.10 获取某天指定时间段详细运动数据

#### 3.11 获取某天的详细睡眠数据

#### 3.12 获取某天指定时间段详细睡眠数据(预留的接口)

#### 3.13 查找戒指(戒指亮灯)

#### 3.14 切换到拍照界面

#### 3.15 保持拍照界面

#### 3.16 停止拍照界面

#### 3.17 硬重启戒指

#### 3.18 获取戒指Mac地址

#### 3.19 获取定时血压测量功能的信息

#### 3.20 设置定时血压测量功能的信息

#### 3.21 获取定时血压测量的历史数据

#### 3.22 重置戒指到出厂设置状态

#### 3.23 获取锻炼历史数据

#### 3.24 获取手动测量血压测量的历史数据

#### 3.25 获取定时心率历史数据

#### 3.26 获取定时心率功能的信息

#### 3.27 设置定时心率功能的信息

#### 3.28 根据指定时间戳, 获取运动数据概要信息

#### 3.29 根据指定新版运动+概要信息, 获取该次运动的部分概要信息和详细数据

#### 3.30 获取/设置用户目标信息

#### 3.31 获取定时体温测量的历史数据

#### 3.32 获取手动体温测量的历史数据

#### 3.33 获取血氧测量的历史数据

#### 3.34 发送固件文件

#### 3.35 收到手环消息

#### 3.36 设置/获取定时血氧开关状态

#### 3.37 发送测量指令 (指令封装在QCSDKManager中)

- 3.38 睡眠协议(获取某一天到今天)
- 3.39 实时心率测量
- 3.40 发起运动
- 3.41 获取压力数据
- 3.42 获取心率变异性(HRV)数据
- 3.43 获取&设置触摸/手势管理
- 3.43 佩戴校准
- 3.44 久坐提醒

## 开发指南

---

### 1 环境配置

#### 1.1 添加framework


将 `QCBandSdk.framework` 加入到项目，framework支持iOS 9.0以上版本

**注意:**因为framework中使用了分类，需要在项目中添加设置

**Target->Build Settings -> Other Linker Flags** 添加 **-ObjC**

在Target-Build Settings-Excluded Architectures中添加以下代码

```
EXCLUDED_ARCHS__EFFECTIVE_PLATFORM_SUFFIX_simulator__NATIVE_ARCH_64_BIT_x86_64=arm64
arm64e armv7 armv7s armv6 armv8 EXCLUDED_ARCHS=$(herited)
$(EXCLUDED_ARCHS__EFFECTIVE_PLATFORM_SUFFIX_$(EFFECTIVE_PLATFORM_SUFFIX)__NATIVE_ARCH_64_
BIT_$(NATIVE_ARCH_64_BIT))
```



## 1.2 配置蓝牙权限

在info.plist文件中配置蓝牙权限

```
<key>NSBluetoothAlwaysUsageDescription</key>
<string>App需要是使用您的蓝牙设备</string>
<key>NSBluetoothPeripheralUsageDescription</key>
<string>App需要是使用您的蓝牙设备</string>
```

## 2 使用

### 2.1 设备支持的服务UUID

在 `QCSDKManager.h` 中定义了，设备支持的服务UUID：

```
extern NSString *const QCBANDSDKSERVERUUID1;
extern NSString *const QCBANDSDKSERVERUUID2;
```

### 2.2 导入framwrok库

在代码中引入framework库

```
#import <QCBandSDK/QCSDKManager.h>
#import <QCBandSDK/QCSDKCmdCreator.h>
```

使用单例初始化 `[QCSDKManager sharedInstance]`

`QCSDKManager`:用于加入连接的外围设备

`QCSDKCmdCreator`:用于向外围设备发送指令

## 2.3 扫描戒指

### 初始化

在权限允许和蓝牙开启状态下，才能开始扫描。

导入Apple的CoreBluetooth库,并遵循两个协议 `<CBCentralManagerDelegate, CBPeripheralDelegate>`

```
#import <CoreBluetooth/CoreBluetooth.h>
```

声明中心角色和外围角色

```
/*中心角色,app*/
@property (strong, nonatomic) CBCentralManager *centerManager;

/*外设角色,扫描到的外围设备*/
@property (strong, nonatomic) NSMutableArray<CBPeripheral *> *peripherals;

/*连接的外设角色*/
@property (strong, nonatomic) CBPeripheral *connectedPeripheral;
```

实例化中心角色

```
self.centerManager = [[CBCentralManager alloc] initWithDelegate:self queue:nil];
```

### 扫描戒指

使用扫描外围设备

```
NSArray *serviceUUIDStrings = @[QCBANDSDKSERVERUUID1,QCBANDSDKSERVERUUID2];

NSMutableArray *uuids = [NSMutableArray array];
for (id obj in serviceUUIDStrings) {
    if ([obj isKindOfClass:[NSString class]]) {
        CBUUID *uuid = [CBUUID UUIDWithString:obj];
        [uuids addObject:uuid];
    }
}

NSDictionary *option = @{CBCentralManagerScanOptionAllowDuplicatesKey : [NSNumber
    numberWithBool:NO]};
[self.centerManager scanForPeripheralsWithServices:uuids options:option];
```

注:在代理中获取扫描到的外围设备,可以通过设备名称等相关信息, 进行二次过滤。

```
- (void)centralManager:(CBCentralManager *)central didDiscoverPeripheral:(CBPeripheral *)peripheral advertisementData:(NSDictionary<NSString *,id> *)advertisementData RSSI:(NSNumber *)RSSI {  
    if (peripheral.name.length > 0) {  
        [self.peripherals addObject:peripheral];  
        [self.deviceList reloadData];  
    }  
}
```

## 2.4 取消扫描设备

调用中心角色的接口, 停止扫描

```
[self.centerManager stopScan];
```

## 2.5 连接戒指

开始连接

```
self.connectedPeripheral = self.peripherals[indexPath.row];  
[self.centerManager connectPeripheral:self.connectedPeripheral options:nil];
```

连接成功后, 传入外围设备到SDK

```
- (void)centralManager:(CBCentralManager *)central didConnectPeripheral:(CBPeripheral *)peripheral {  
    [[QCSDKManager sharedInstance] addPeripheral:peripheral];  
}
```

## 2.6 断开连接

```
[self.centerManager cancelPeripheralConnection:self.connectedPeripheral];
```

断开连接后, 删除外围设备

```

- (void)centralManager:(CBCentralManager *)central didDisconnectPeripheral:(CBPeripheral
*)peripheral error:(nullable NSError *)error {
    [[QCSDKManager sharedInstance] removePeripheral:peripheral];
}

```

## 3. 戒指支持的指令

### 3.1 设置手环时间

```

/**
 * Set the time of the watch
 * 设置戒指的时间
 *
 */
+ (void)setTime:(NSDate *)date success:(void (^)(NSDictionary *featureList))suc failed:
(void (^)(void))fail;

```

### 3.2 读取戒指电量

```

/*!
 * @func 读取戒指电量
 * @param suc battery:电量等级(0~8)
 */
+ (void)readBatterySuccess:(void (^)(int battery))suc failed:(void (^)(void))fail;

```

### 3.3 绑定戒指亮灯

```

/**
 * 绑定戒指亮灯
 */
+ (void>alertBindingSuccess:(nullable void (^)(void))suc fail:(nullable void (^)(
void))fail;

```

### 3.4 设置戒指时间进制/用户个人信息

```

/**
 设置戒指时间进制/用户个人信息
 @param twentyfourHourFormat      : YES 24小时制; NO 12小时制
 @param metricSystem              : YES 公制; NO 英制
 @param gender                    : 性别 (0=男, 1=女)
 @param age                       : 年龄 (岁)

```

```

@param height          : 身高 (厘米)
@param weight          : 体重 (kg)
@param sbpBase         : 收缩压基值 (mmhg) (保留值, 可设置0)
@param dbpBase         : 舒张压基值 (mmhg) (保留值, 可设置0)
@param hrAlarmValue    : 心率报警值 (bpm) (保留值, 可设置0)
*/
+ (void)setTimeFormatTwentyfourHourFormat:(BOOL)twentyfourHourFormat
    metricSystem:(BOOL)metricSystem
    gender:(NSInteger)gender
    age:(NSInteger)age
    height:(NSInteger)height
    weight:(NSInteger)weight
    sbpBase:(NSInteger)sbpBase
    dbpBase:(NSInteger)dbpBase
    hrAlarmValue:(NSInteger)hrAlarmValue
    success:(void (^)(BOOL, BOOL, NSInteger, NSInteger, NSInteger, NSInteger, NSInteger,
NSInteger, NSInteger))success
    fail:(void (^)(void))fail;

```

### 3.5 获取戒指时间进制/用户个人信息

```

/**
获取戒指时间进制/用户个人信息
@param success isTwentyfour: YES 24小时制; NO 12小时制
            isMetricSystem: YES 公制; NO 英制
            gender: 性别 (0=男, 1=女)
            age: 年龄 (岁)
            height: 身高 (厘米)
            weight: 体重 (kg)
            sbpBase: 收缩压基值 (mmhg) (保留值, 可设置0)
            dbpBase: 舒张压基值 (mmhg) (保留值, 可设置0)
            hrAlarmValue: 心率报警值 (bpm) (保留值, 可设置0)
*/
+ (void)getTimeFormatInfo:(nullable void (^)(BOOL isTwentyfour, BOOL isMetricSystem,
NSInteger gender, NSInteger age, NSInteger height, NSInteger weight, NSInteger sbpBase,
NSInteger dbpBase, NSInteger hrAlarmValue))success fail:(nullable void (^)(void))fail;

```

### 3.6 获取戒指固件的版本号

```

/**
* @func 获取戒指固件(Application)的版本号
* @param success 软件硬件版本号格式一般为"x.x.x"
*/
+ (void)getDeviceSoftAndHardVersionSuccess:(void (^)(NSString *_Nonnull, NSString
*_Nonnull))success fail:(void (^)(void))fail;

```

### 3.7 获取当前计步信息

```
/*!
 * @func    获取当前计步信息(可以同步最新纪录)
 */
+ (void)getCurrentSportSucess:(void (^)(SportModel *sport))suc failed:(void (^)(void))fail;
```

### 3.8 获取某天总的统计数据

```
/*!
 * @func    获取某天总的统计数据
 * @param   index    : 0->当天 1->1天前.最大29
 * @param   suc      : 使用这个命令不能准确的获取当天的统计数据,设备没15分钟会存一次数据,所以会有15分钟间隔
 */
+ (void)getOneDaySportBy:(NSInteger)index success:(void (^)(SportModel *model))suc fail:(void (^)(void))fail;
```

### 3.9 获取某天的详细运动数据

```
/*!
 * @func    获取某天的详细运动数据
 * @discussion 每15分钟一个刻度, 每天最多会有96条数据. 详细请看返回内容
 * @param   items    sports:返回所有的运动model
 */
+ (void)getSportDetailDataByDay:(NSInteger)dayIndex sportDatas:(nullable void (^)(NSArray<SportModel *> *sports))items fail:(nullable void (^)(void))fail;
```

### 3.10 获取某天指定时间段详细运动数据



```

/*!
 * @func 获取某天指定时间段详细运动数据
 * @param minuteInterval 每个索引的分钟间隔
 * @param beginIndex 时间段开始索引
 * @param endIndex 时间段结束索引
 * @param items sports:返回所有的运动model
 */
+ (void)getSportDetailDataByDay:(NSInteger)dayIndex minuteInterval:
(NSInteger)minuteInterval beginIndex:(NSInteger)beginIndex endIndex:(NSInteger)endIndex
sportDatas:(nullable void (^)(NSArray<SportModel *> *sports))items fail:(nullable void
(^)(void))fail;

```

### 3.11 获取某天的详细睡眠数据

```

//睡眠数据类型
typedef NS_ENUM(NSInteger, SLEEPTYPE) {
    SLEEPTYPENONE,    //无数据
    SLEEPTYPESOBER,   //清醒
    SLEEPTYPELIGHT,   //浅睡
    SLEEPTYPEDEEP,    //深睡
    SLEEPTYPEUNWEARED //未佩戴
};

@interface QCSleepModel : NSObject
@property (nonatomic, assign) SLEEPTYPE type;           //睡眠类型
@property (nonatomic, strong) NSString *happenDate;    //发生时间 yyyy-MM-dd HH:mm:ss
@property (nonatomic, strong) NSString *endTime;       //结束时间 yyyy-MM-dd HH:mm:ss.
@property (nonatomic, assign) NSInteger total;         //开始时间与结束时间的时间间隔(单位: 分钟)
@end

/*!
 * @func 获取某天的详细睡眠数据
 * @discussion 每种睡眠类型对应的时间段, 详细请看返回内容
 * @param items sleeps:返回所有的睡眠model
 */
+ (void)getSleepDetailDataByDay:(NSInteger)dayIndex sleepDatas:(nullable void (^)(
NSArray<QCSleepModel *> *sleeps))items fail:(nullable void (^)(void))fail;

```

### 3.12 获取某天指定时间段详细睡眠数据(预留的接口)

```

/*!
 * @func    获取某天指定时间段详细睡眠数据(预留的接口)
 * @param  minuteInterval  每个索引的分钟间隔
 * @param  beginIndex      时间段开始索引
 * @param  endIndex        时间段结束索引
 * @param  items            sports:返回所有的睡眠model
 */
+ (void)getSleepDetailDataByDay:(NSInteger)dayIndex minuteInterval:
(NSInteger)minuteInterval beginIndex:(NSInteger)beginIndex endIndex:(NSInteger)endIndex
sleepDatas:(nullable void (^)(NSArray<SleepModel *> *sleeps))items fail:(nullable void
(^)(void))fail;

```

### 3.13 查找戒指(戒指亮灯)

```

/**
 * 查找戒指
 */
+ (void)lookupDeviceSuccess:(void (^)(void))suc fail:(void (^)(void))fail;

```

### 3.14 切换到拍照界面

```

/**
 * 切换下位机(手环)到拍照界面
 */
+ (void)switchToPhotoUISuccess:(nullable void (^)(void))success fail:(nullable void (^)(void))fail;

```

### 3.15 保持拍照界面

```

/**
 * 保持下位机(手环)拍照界面
 */
+ (void)holdPhotoUISuccess:(nullable void (^)(void))success fail:(nullable void (^)(void))fail;

```

### 3.16 停止拍照界面

```

/**
 * 停止下位机(手环)拍照
 */
+ (void)stopTakingPhotoSuccess:(nullable void (^)(void))success fail:(nullable void (^)(void))fail;

```

### 3.17 硬重启戒指

```

/**
 硬重启戒指
 */
+ (void)resetBandHardlySuccess:(nullable void (^)(void))suc fail:(nullable void (^)(void))fail;

```

### 3.18 获取戒指Mac地址

```

/**
 * @func 获取戒指Mac地址
 * @param success Mac地址格式为"AA:BB:CC:DD:EE:FF"
 */
+ (void)getDeviceMacAddressSuccess:(nullable void (^)(NSString *_Nullable macAddress))success fail:(nullable void (^)(void))fail;

```

### 3.19 获取定时血压测量功能的信息

```

/**
 * 获取定时血压测量功能的信息
 * @param success featureOn YES: 开启; NO: 关闭
 *                  beginTime 开始时间, 格式为"HH:mm"
 *                  endTime 结束时间, 格式为"HH:mm"
 *                  minuteInterval 分钟间隔
 */
+ (void)getSchedualBPInfo:(nullable void (^)(BOOL featureOn, NSString *beginTime, NSString *endTime, NSInteger minuteInterval))success fail:(void (^)(void))fail;

```

### 3.20 设置定时血压测量功能的信息

```

/**
 * 设置定时血压测量功能的信息
 * @param featureOn YES: 开启; NO: 关闭
 * @param beginTime 开始时间, 格式为"HH:mm"
 * @param endTime 结束时间, 格式为"HH:mm"
 * @param minuteInterval 分钟间隔
 */
+ (void)setSchedualBPInfoOn:(BOOL)featureOn beginTime:(NSString *)beginTime endTime:
(NSString *)endTime minuteInterval:(NSInteger)minuteInterval success:(nullable void (^)(
(BOOL featureOn, NSString *beginTime, NSString *endTime, NSInteger
minuteInterval)))success fail:(void (^)(void))fail;

```

### 3.21 获取定时血压测量的历史数据

```

/**
 * 获取定时血压测量的历史数据
 * @param userAge :用户年龄
 * @param success data 心率模块数据, 目前的回复其实统一为心率, 回调中可自行处理
 */
+ (void)getSchedualBPHistoryDataWithUserAge:(NSInteger)userAge success:(nullable void (^)(
NSArray<BloodPressureModel *> *data))success fail:(nullable void (^)(void))fail;

```

### 3.22 重置戒指到出厂设置状态

```

/**
 * 重置戒指到出厂设置状态, 慎用
 */
+ (void)resetBandToFacotrySuccess:(nullable void (^)(void))success fail:(nullable void
(^)(void))fail;

```

### 3.23 获取锻炼历史数据

```

/**
 * @func 获取锻炼历史数据
 * @param lastUnixSeconds 最后一条锻炼数据的发生时间(距1970-01-01 00:00:00的秒数)
 * @note success models 锻炼数据数组
 */
+ (void)getExerciseDataWithLastUnixSeconds:(NSUInteger)lastUnixSeconds getData:(nullable
void (^)(NSArray<ExerciseModel *> *models))getData fail:(nullable void (^)(void))fail;

```

### 3.24 获取手动测量血压测量的历史数据

```
/**
 * 获取手动测量血压测量的历史数据
 * @param lastUnixSeconds 最后一条手动血压数据的发生时间(距1970-01-01 00:00:00的秒数)
 * @param success data 血压数据数组
 */
+ (void)getManualBloodPressureDataWithLastUnixSeconds:(NSUInteger)lastUnixSeconds
success:(nullable void (^)(NSArray<BloodPressureModel *> *data))success fail:(nullable
void (^)(void))fail;
```

### 3.25 获取定时心率历史数据

```
/**
 * @func 获取定时心率历史数据
 * @param dates 需要获取历史数据的日期列表
 * @note success models 定时心率数据数组
 */
+ (void)getSchedualHeartRateDataWithDates:(NSArray<NSDate *> *)dates success:(nullable
void (^)(NSArray<SchedualHeartRateModel *> *models))success fail:(nullable void (^)(
void))fail;

/**
 * @func 获取定时心率历史数据
 * @param dayIndexs 需要获取历史数据的天数(0->今天, 1->昨天, 2->前天, 依次类推)
 * @note success models 定时心率数据数组
 */
+ (void)getSchedualHeartRateDataWithDayIndexs:(NSArray<NSNumber*> *)dayIndexs success:
(void (^)(NSArray<QCSchedualHeartRateModel *> *_Nonnull))success fail:(void (^)(
void))fail;
```

### 3.26 获取定时心率功能的信息

```
/**
 * 获取定时心率功能的信息
 * @param success enable 定时心率功能是否开启. YES: 开启; NO: 关闭
 */
+ (void)getSchedualHeartRateStatusWithCurrentState:(BOOL)enable success:(nullable void
(^)(BOOL enable))success fail:(nullable void (^)(void))fail;
```

### 3.27 设置定时心率功能的信息

```
/**
 * 设置定时心率功能的信息
 * @param enable 定时心率功能是否开启。YES：开启；NO：关闭
 */
+ (void)setSchedualHeartRateStatus:(BOOL)enable success:(nullable void (^)(BOOL enable))success fail:(nullable void (^)(void))fail;
```

### 3.28 根据指定时间戳, 获取运动数据概要信息

```
/**
 根据指定时间戳, 获取该时间戳后的新版运动+(v2)数据概要信息
 @param timestamp 时间戳
 @param finished spSummary - 运动+概要信息数组
 */
+ (void)getSportPlusSummaryFromTimestamp:(NSTimeInterval)timestamp finished:(nullable void (^)(NSArray *_Nullable spSummary, NSError *_Nullable error))finished;
```

### 3.29 根据指定新版运动+概要信息, 获取该次运动的部分概要信息和详细数据

```
/**
 根据指定新版运动+概要信息, 获取该次运动的部分概要信息和详细数据
 @param finished spSummary - 运动+概要信息数组
 */
+ (void)getSportPlusDetailsWithSummary:(OdmGeneralExerciseSummaryModel *)summary
finished:(nullable void (^)(OdmGeneralExerciseSummaryModel *_Nullable summary,
OdmGeneralExerciseDetailModel *_Nullable detail, NSError *_Nullable error))finished;
```

### 3.30 获取/设置用户目标信息

```
/**
 获取用户目标信息
 @param suc stepTarget 步数目标
        calorieTarget 卡路里目标, 单位: 卡
        distanceTarget 距离目标, 单位: 米
        sportDuration 运动时长目标 单位: 分钟(保留值, 默认: 0)
        sleepDuration 睡眠时长目标 单位: 分钟(保留值, 默认: 0)
 */
```

```

+ (void)getStepTargetInfoWithSuccess:(nullable void (^)(NSInteger stepTarget,NSInteger
calorieTarget,NSInteger distanceTarget,NSInteger sportDuration,NSInteger
sleepDuration))suc fail:(nullable void (^)(void))fail;

/**
 设置用户目标信息
  @param stepTarget 步数目标
  @param calorieTarget 卡路里目标, 单位: 卡
  @param distanceTarget 距离目标, 单位: 米
  @param sportDuration 运动时长目标 单位: 分钟(保留值, 默认: 0)
  @param sleepDuration 睡眠时长目标 单位: 分钟(保留值, 默认: 0)
  */
+ (void)setStepTarget:(NSInteger)stepTarget calorieTarget:(NSInteger)calorieTarget
distanceTarget:(NSInteger)distanceTarget sportDurationTarget:(NSInteger)sportDuration
sleepDurationTarget:(NSInteger)sleepDuration success:(nullable void (^)(void))suc fail:
(nullable void (^)(void))fail;

```

### 3.31 获取定时体温测量的历史数据

```

/**
 * 获取定时体温测量的历史数据
 */
+ (void)getSchedualTemperatureDataByDayIndex:(NSInteger)dayIndex finished:(nullable void
(^)(NSArray *_Nullable temperatureList, NSError *_Nullable error))finished;

```

### 3.32 获取手动体温测量的历史数据

```

/**
 * 获取手动体温测量的历史数据
 */
+ (void)getManualTemperatureDataByDayIndex:(NSInteger)dayIndex finished:(nullable void
(^)(NSArray *_Nullable temperatureList, NSError *_Nullable error))finished;

```

### 3.33 获取血氧测量的历史数据

```

/**
 * 获取血氧测量的历史数据
 */
+ (void)getBloodOxygenDataByDayIndex:(NSInteger)dayIndex finished:(void (^)(NSArray *
_Nullable, NSError *_Nullable))finished;

```

### 3.34 发送固件文件

```
/**
 发送固件文件，要求使用bin文件升级，结果会在回调里边处理

  @param data          OTA二进制字符流
  @param start          开始发送回调
  @param percentage     进度回调
  @param success        成功回调
  @param failed         失败回调
  */

+ (void)syncOtaBinData:(NSData *)data
    start:(nullable void (^)(void))start
  percentage:(nullable void (^)(int percentage))percentage
    success:(nullable void (^)(int seconds))success
    failed:(nullable void (^)(NSError *error))failed;
```

### 3.35 收到手环消息

```
@interface QCSDKManager : NSObject

/**
 * 收到手表，查找手机的通知
 */
@property(nonatomic,copy)void(^findPhone)(void);

/**
 * 收到手表，进入相机的通知
 */
@property(nonatomic,copy)void(^switchToPicture)(void);

/**
 * 收到手表，拍照的通知的回调
 */
@property(nonatomic,copy)void(^takePicture)(void);

/**
 * 收到手表，结束拍照的通知的回调
 */
@property(nonatomic,copy)void(^stopTakePicture)(void);

// 单例类实例
+ (instancetype)shareInstance;

@end
```



### 3.36 设置/获取定时血氧开关状态

```
/**
 * 设置定时血氧测量功能的信息
 * @param featureOn YES: 开启; NO: 关闭
 */

+ (void)setSchedualBOInfoOn:(BOOL)featureOn success:(nullable void (^)(BOOL featureOn))success fail:(void (^)(void))fail;

/**
 * 获取定时血氧测量功能的信息
 * @param success featureOn YES: 开启; NO: 关闭
 */

+ (void)getSchedualBOInfoSuccess:(nullable void (^)(BOOL featureOn))success fail:(void (^)(void))fail;
```

### 3.37 发送测量指令 (指令封装在QCSDKManager中)

```
typedef NS_ENUM(NSInteger, QCMeasuringType) {
    QCMeasuringTypeHeartRate = 0, //Heart rate measurement
    QCMeasuringTypeBloodPressue, //blood pressure measurement
    QCMeasuringTypeBloodOxygen, //blood oxygen measurement
    QCMeasuringTypeOneKeyMeasure, //One-click measurement
    QCMeasuringTypeStress,
    QCMeasuringTypeBloodGlucose,
    QCMeasuringTypeCount,
};

//测量结果为hanle回调中的result
//测量心率的时候,result返回的是NSNumber: @(60)
//测量血压的时候,result返回的是NSDictionary:@{@"sbp":@"120",@"dbp":@"60"}
//测量血氧的时候,result返回的是NSNumber: @(98)

/// Send measurement order
/// @param type :Measurement type
/// @param measuring :Real-Time Measuring Value
/// @param handle :Measurement result callback (error code: -1: failed to send start command, -2: failed to send end command, -3: bracelet is not properly worn)
- (void)startToMeasuringWithOperateType:(QCMeasuringType)type measuringHandle:(void (^)(id _Nullable result))measuring completedHandle:(void (^)(BOOL isSuccess,id _Nullable result, NSError * _Nullable error))handle;

/// Send measurement order
/// @param type :Measurement type
```

```

/// @param measuring      :Real-Time Measuring Value
/// @param handle         :Measurement result callback (error code: -1: failed to
send start command, -2: failed to send end command, -3: bracelet is not properly worn)
- (void)startToMeasuringWithOperateType:(QCMeasuringType)type timeout:(NSInteger)timeout
measuringHandle:(void (^)(id _Nullable result))measuring completedHandle:(void (^)(BOOL
isSuccess,id _Nullable result, NSError * _Nullable error))handle;

/// stop measurement command
/// @param type           :Measurement type
/// @param handle         :Measurement result callback (error code:-1: Failed to send end
command)
- (void)stopToMeasuringWithOperateType:(QCMeasuringType)type completedHandle:(void (^)(
BOOL isSuccess, NSError *error))handle;

```

### 3.38 睡眠协议(获取某一天到今天)

```

/*!
 * @func    获取从某天到今天的所有睡眠数据
 * @param fromDayIndex  距离今天的天数, (0:表示今天, 1: 表示昨天)
 * @param items        返回的睡眠数据(key: 距离今天的天数, value: 对应的睡眠数据)
 * @param fail         失败的 回调
 */
+ (void)getSleepDetailDataFromDay:(NSInteger)fromDayIndex sleepDatas:(nullable void (^)(
NSDictionary <NSString*, NSArray<QCSleepModel*>*> _Nonnull))items fail:(nullable void (^)(
void))fail;

```

### 3.39 实时心率测量

```

typedef enum {
    QCBandRealTimeHeartRateCmdTypeStart = 0x01, //Start real-time heart rate measurement
    QCBandRealTimeHeartRateCmdTypeEnd, //End real-time heart rate measurement
    QCBandRealTimeHeartRateCmdTypeHold, //Continuous heart rate test (for continuous
measurement to keep alive)
} QCBandRealTimeHeartRateCmdType;

/**
 * RealTime HeartRate Measuring
 * 实时心率测量
 *
 * @param type      :commond type
 * @param finished  :finish callback
 */
+ (void)realTimeHeartRateWithCmd:(QCBandRealTimeHeartRateCmdType)type finished:(nullable
void (^)(BOOL))finished;

```

### 3.40 发起运动

```
/// Set Sport Mode State
///
/// - Parameters:
///   - sportType: type
///   - state: state
///   - finished: finished callback
+ (void)operateSportModeWithType:(OdmSportPlusExerciseModelType)sportType state:
(QCSportState)state finish:(void (^)(id _Nullable, NSError * _Nullable))finished;
```

获取回调:

```
[QCSDKManager sharedInstance].currentSportInfo = ^(QCSportInfoModel * _Nonnull
sportInfo) {

    NSLog(@"sportType:%zd,duration:%zd,state:%u,hr:%zd,step:%zd,calorie(unit:calorie):%zd,di
stance(unit:meter):%zd",sportInfo.sportType,sportInfo.duration,sportInfo.state,sportInfo.
hr,sportInfo.step,sportInfo.calorie,sportInfo.distance);
};
```

### 3.41 获取压力数据

```
/// Get Schedual Stress Datas (Only Ring Support)
///
/// - Parameters:
///   - dates: 0-6,0:today,1:yesterday....
///   - finished: finished callback
+ (void)getSchedualStressDataWithDates:(NSArray<NSNumber*> *)dates finished:(void (^)(
NSArray * _Nullable, NSError * _Nullable))finished;;

/// Get Schedual Stress Status
///
/// - Parameter finished: finished callback
+ (void)getSchedualStressStatusWithFinshed:(nullable void (^)(BOOL,NSError * _Nullable
error))finished;

/// Set Schedual Stress Status
///
/// - Parameters:
///   - enable:YES:On,NO:Off
///   - finished: finished callback
+ (void)setSchedualStressStatus:(BOOL)enable finshed:(nullable void (^)(NSError
*_Nullable error))finished;
```

### 3.42 获取心率变异性(HRV)数据

```
/// Get Scheduling HRV Data (Only Ring Support)
///
/// @param dates 0-6,0:today,1:yesterday....
/// @param finished finished callback
+ (void)getSchedulingHRVDataWithDates:(NSArray<NSNumber*> *)dates finished:(void (^)(
    NSArray * _Nullable, NSError * _Nullable))finished;

/// Get Scheduling HRV Status
///
/// - Parameter finished: finished callback
+ (void)getSchedulingHRVWithFinished:(nullable void (^)(BOOL, NSError * _Nullable
    error))finished;

/// Set Scheduling HRV Status
///
/// - Parameters:
///   - enable:YES:On,NO:Off
///   - finished: finished callback
+ (void)setSchedulingHRVStatus:(BOOL)enable finished:(nullable void (^)(NSError * _Nullable
    error))finished;
```

### 3.43 获取&设置触摸/手势管理

```
/// Get Touch Control Type
///
/// @param finished : callback-> type:QCTouchGestureControlType ,strength:1-10
+ (void)getTouchControlFinished:(nullable void (^)(
    QCTouchGestureControlType, NSInteger, NSError * _Nullable error))finished;

/// Set Touch Control Type
/// @param type : type
/// @param strength :1-10
/// @param finished :callback
+ (void)setTouchControl:(QCTouchGestureControlType)type strength:(NSInteger)strength
    finished:(nullable void (^)(NSError * _Nullable error))finished;

/// Get Gesture Control Type
///
/// @param finished : callback-> type:QCTouchGestureControlType ,strength:1-10
+ (void)getGestureControlFinished:(nullable void (^)(
    QCTouchGestureControlType, NSInteger, NSError * _Nullable error))finished;

/// Set Gesture Control Type
```

```

/// @param type : type
/// @param strength :1-10
/// @param finished :callback
+ (void)setGestureControl:(QCTouchGestureControlType)type strength:(NSInteger)strength
finished:(nullable void (^)(NSError *_Nullable error))finished;

```

### 3.43 佩戴校准

```

/// Wearing Calibration
///
/// @param type 1->Start calibration (reset ring data), 2->End calibration, 3->Get single
data, 4->Power consumption mode, 5->Stop power consumption, 6->App starts calibration
/// @param finished finished callback
+ (void)wearCalibration:(NSInteger)type finished:(nullable void (^)(NSError *_Nullable
error))finished;

```

### 3.44 久坐提醒

```

/// Get Sedentary Reminder (Only Ring Support)
///
/// @param fromDayIndex :0->Today,1->Yesterday,2->The day before yesterday ....
/// @param finished : callback
+ (void)getSedentaryReminderFromDay:(NSInteger)fromDayIndex finished:(nullable void (^)(
NSDictionary <NSString*,NSArray<QCSedentaryModel*>*>*_Nullable datas, NSError *_Nullable
error))finished;

```